


## Flow Cytometry Adaptor & Calibration Slide



Cell Sorter Adaptor with AmpliGrid AG480F



Cell Sorter Adaptor with Calibration Slide



AmpliSpeed slide cycler

### Features

- High precision SBS size adaptor
- Accurate placement of AmpliGrid slides
- Fits in various commercial flow cytometry sorters
- Positioning accuracy +/- 50  $\mu$ m

### Features Calibration Slide

- Metal plated AmpliGrid slide
- Optimized visualization of all 48 reaction sites
- Ideal tool for setting sorting positions
- Increased sorting accuracy

### Advantages

- Flow Cytometry Sorter setup for AmpliGrid made easy
- Quick starting guides available for
- Beckman Coulter MoFlo™
- BD FACSVantage™
- BD FACSAria™

### Order Information

Ref.-No	Product	Description
OAX04118	Flow Cytometry Adaptor	Adaptor for precise placement of one AG480F
OAX04120	Flow Cytometry Calibration Slide	AG480F metalized for Flow Cytometry Sorter calibration

BD, the BD log, FACS™, FACSAria™ and FACSFlow™ are trademarks of Becton, Dickinson and Company.

For research use only. Not for use in diagnostic procedures.

[www.advalytix.com](http://www.advalytix.com)