

Microarray coverslips with precision spacer



AdvaGap AP10X

Features

- precision spacer of 50 μm thickness
- different sizes
- controlled glass quality
- consistent flatness
- low parasitic DNA-binding

Applications

- oligo microarrays
- cDNA microarrays
- protein microarrays

Benefits

- uniform spreading of hybridization solution
- no air bubbles
- low background
- homogeneous hybridization

Order Information

Ref.-No	Description
OAX05212	AP101-50: cover size 22x22 mm with approx. volume 22 μL (50 Pcs/pack)
OAX05213	AP102-25: cover size 22x40 mm with approx. volume 40 μL (25 Pcs/pack)
OAX05214	AP103-25: cover size 22x50 mm with approx. volume 50 μL (25 Pcs/pack)
OAX05215	AP104-25: cover size 22x60 mm with approx. volume 60 μL (25 Pcs/pack)

For research use only. Not for use in diagnostic procedures.