



AdvaGold & AdvaBlue PCR dyes



AdvaBlue / AdvaGold



AmpliGrid AG480F



No influence on C_t values

Features

- Additive dye for PCR / RT-PCR master mixes
- Colors aqueous phase of assay for easy distinction from sealing solution
- Choice of blue or gold dyes to avoid conflict with fluorescent detection
- Biologically inert

Applications

- Use to more easily distinguish the aqueous phase (PCR product) from sealing solution when retrieving the amplicon post PCR
- Compatible with all AmpliGrid applications (standard and multiplex RT-PCR, PCR, sequencing and qPCR pre-amplification)

Advantages

- Improved visualization of small volumes
- Helps avoid handling errors during pipetting assay and sealing solution on and off AmpliGrid slides
- No influence on C_t values during qPCR

Order Information

Ref.-No	Product	Description
OAX04243	AdvaBlue	10x concentrated PCR dye
OAX04227	AdvaGold	10x concentrated PCR dye

For research use only. Not for use in diagnostic procedures.

www.advalytix.com/RT-PCR