

Necessary consumables for the use of the SlideBooster

AdvaHum AM10X

Features

- guarantees 100% humidity with most commercially sold or selfmade hybridization buffers
- low evaporation of hybridization buffer

Benefits

- high signal to noise ratio



AdvaHum AM101 and AM102


AdvaSon AS100

Features

- builds a uniform layer of liquid underneath the slide
- transfers SAW from chip to slide

Benefits

- homogeneous agitation



AdvaSon AS100

Order Information

Ref.-No	Description
OAX05210	AM101: humidifying solution for formamid-based hybridization buffer (1 x 50 mL)
OAX05211	AM102: humidifying solution for formamid-free hybridization buffer (1 x 50 mL)
OAX05209	AS100: coupling liquid (1 x 50 mL)

For research use only. Not for use in diagnostic procedures.