


Slide processing tube for the AdvaWash station



AdvaWash AW400 (Advalytix) including Advatubes



Placing the AdvaTube in the AdvaWash chamber



Slide positions in the AdvaTube

Features

- up to 4 slides sized 24.3 - 26 mm x 75-76 mm
- volume 50 mL
- compatible with centrifuge holders for „Falcon“ tubes
- reusable for up to 10 runs

Applications

- microarrays: oligo, cDNA and protein
- CGH: oligo and BAC
- immunohistochemistry procedures
- Slide-based FISH procedures

Benefits

- optimized geometry for maximum washing efficiency
- special fluidics for flawless drying
- tinted to protect sensitive dye markers from bleaching
- high chemical resistivity
- reduction of manual handling steps

Order Information

Ref.-No	Description
OAX05216	AdvaTube: slide processing tube (10 pcs/pack)
OAX05217	AdvaTube: slide processing tube (50 pcs/pack)
OAX05218	AdvaTube: slide processing tube (100 pcs/pack)

For research use only. Not for use in diagnostic procedures.