


## SAW sample mixer for BD FACSAria™



CellEvatorAria sample mixer module

### Features

- Powered by Advalytix surface acoustic wave (SAW) mixing technology
- Mixes sample from the bottom up to prevent cell settlement
- Gentle mixing technology - no adverse effects on cells
- Easy installation in the BD FACSAria™
- Compatible with 5 ml and 15 ml BD polystyrene sample tubes


### Applications

Prevents drop-off of flow rates and reduces clogging for:


- long sorts
- sorting large or fragile cells
- sorting primary or immortalized cells
- sorting bacterial suspensions

### Advantages

- Prevents flow rate drop-off and clogging more effectively than built-in mixer
- Constant flow rate up to the last sample droplet, higher yield
- Vertical mixing from the bottom up most effective to prevent cell settlement
- Reduced cell necrosis



CellEvator mixer module installed in BD FACSAria™



Comparison of sorting performance

### Order Information

Ref.-No.	Product	Description
OAX06160	CellEvatorAria	CellEvatorAria, acoustic sample agitator for BD FACSAria™ cell sorters

BD, the BD logo, FACS™, FACSAria™ and FACSFlow™ are registered trademarks of Becton, Dickinson and Company.